
*Bob
AND
Brad®*

User Manual

QL/MINI.Q2-D
Mini Massage Gun

Warnings and Safety Cautions

It is important to know the following before use, otherwise it may cause damage or become dangerous!

- For use only by healthy adults other than older adults;
- Use the product only on the soft tissues of the body without causing pain or discomfort;
- **Do not** use on the head, or any hard or bony part of the body;
- Either controlled downward pressure or the application of pressure may produce bruising. Check the massage area frequently and - stop using the product as soon as signs of pain or discomfort appear;
- Use only on dry, clean skin surfaces on body parts, pressing and moving gently on the skin;
- Keep your fingers, hair or other body parts away from the shaft and back of the vibrating head, as they may get caught;
- **Do not** apply to bony areas where there is less muscle coverage, such as the shins, backs of the feet, back of the hands;
- **Do not** use the product on or around sensitive areas such as eyes, teeth, breast implants, genital organs, artificial limbs, etc;
- **Do not** place the product in any area where there is a risk of fracture;
- If the pain is persistent, discontinue use immediately and consult a physician.
- This product is not a medical device, it's not

- recommended for medical use, not intended for medical use in hospitals;
- Please check the main unit and cable before using, if damaged, it must be replaced by the manufacturer or its appointed professional to avoid danger;

FCC Statement

1.This device complies with Part 15 of the FCC Rules. Operation is subject to the following.

Two conditions:

- (1) This device may not cause harmful interference.
- (2) This device must accept any interference received, including interference that may cause undesired operation.

FCC Statement NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct

the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help with the supplied fastener. **DO NOT** hang the unit by the power cord.

Disclaimer


Please contact your physician before using the **Bob and Brad QL/MINI.Q2-D massage gun** if you have any physical health concerns or have any of the following conditions:

Pregnancy, diabetes, epilepsy, migraines, herniated disks, spondylolisthesis, recent joint replacements or IUDs. Please also contact your physician before use if you have undergone recent surgery or have any metal pins or plates in your body.

Contact Us

www.bobandbrad.com
support@bobandbrad.com
Tel:+1 (612)567-3035
M-F 9:00AM-5:30PM

Main Structure and The Fittings

Main Structure and The Fittings

1. Fork Head

2. Ball Head

3. Bullet Head

4. Air Cushion Head


5. Flat Head

Operation Steps

Removing And Replacing A Massage Head

1. Before attaching the massage head, turn off the massage gun.
2. Grab the inserted head and firmly pull out to remove it.
3. Insert the stem of a new head into the unit, press firmly to secure.

(Note: Fork-head need to be installed as shown in the attached picture, otherwise there is a risk of falling off!)

Turn the device ON/OFF And Changing Percussion Massage Frequency

1. Turn the device ON/OFF In the shutdown state, press the "ⓘ" button on the panel for 2 seconds to turn on

2. To changing percussion massage frequency when the device is on, press the button "ⓘ" to turn on the first gear. If you would like to advance to a higher setting, press the button "ⓘ" again. The active gear will be indicated by the back light. To pause the massage unit the button "ⓘ" for 2 seconds. Press for 2 seconds again to resume.

Noted:

- A. After 10 minutes of operation, the device automatically changes to level 0 and stops output;
- B. The LED light indicates the corresponding strike frequency level, the lowest at level 1 and the highest at level 5;
- C. When the device is in standby (0 level), it will automatically shut down without any operation after 20 seconds.

Battery Display

When the device is running, the battery indicator is red, indicating that the power is less than 30%; when the indicator is red and blinking it's indicating that the power is less than 20%, and should be charged soon.

Charging

1. Before the first use, please fully charge the device. To charge fully, it takes 3.5 hours.
2. Connect the output end of the connection line to a DC charging port at the bottom of the device, and then connect the input end of the connection

line to the mobile power source or mobile phone charger (5V/2A, if the charging current or voltage is less than 5 V/2 A, charging time will be longer. Any regular phone adapters will be able to charge the massage gun without a problem.)

3. When connecting the power cable the light will appear green and, it means the device is charging, when fully charged the green indicator will turn off.
4. The battery can be charged at any power level; when charging, please ensure that the device is off, when the device is running, charging will automatically be protected.

Technical Details:

QL/MINI.Q2-D Mini Massage Gun

Speed: 1800-3000 r/min

Rated Power: 25W

Dimensions: 14.6x8.6x4.8cm / 5.7"x3.3"x1.8"

Battery:

Rated Voltage: 11.1V 

Battery capacity: 2500mAh

Charging time: 210 minutes

Guarantee

This product comes with a one-year warranty from the date of purchase. Warranty will be void due to any of the following:

Damaged due to improper usage or storage.

Damaged by attempting to repair or disassembling the device. Period of warranty exceeded.

Utilize the Bob and Brad massage gun before, during, or after your workout for maximum benefits. This massage gun can help accelerate muscle tissue responsiveness by increasing hydration, range of motion and function to tight muscles.

Indication	Action	Speed	Pressure	Indication	Duration	Muscle Position	Movement
Warm Up / Activation	Activate and prime muscles prior to exercise. Reduce risk of injury.	Level 3	Light to Moderate	Warm Up /Activation	30 sec or less	Relaxed. Move along the muscle.	Quick, high speed over area.
Mid Workout / Myofascial Mobility	Help muscle recovery and reactivation.	Level 4	Heavy (especially for large muscle groups)	Mid Workout / Myofascial Mobility	8-12 sec	Target specific area.	Localised
Recovery	Relax muscles and aid recovery. Down-regulate or turn off muscles.	Level 1 or 2	Moderate to Heavy	Recovery	90 sec or more	Move along the muscle.	Slow movement along the muscle
Pain Modulation	Reduction in pain to increase movement.	Level 1	Light	Pain Modulation	10 to 30 sec	Relaxed. Localised spot.	Over painful area.
Body Awareness	Improve movement or mechanics of a joint. Activation of a muscle. Draw attention to a muscle	All Levels	Light to Moderate	Body	10 to 20 sec	Relaxed or with movement	Over specific target area or movement over area.

EN

DE


FR

ES

IT

Full Body Massage



NECK

The ball, air cushion and flat head attachments are recommended for neck and shoulder muscles. Please avoid using on the skull. We recommend to use the air cushion attachment for thin or less muscled areas to avoid pain or discomfort.

Full Body Massage

ARM & PLAM

The ball, air cushion and flat head attachments are recommended for massaging the arm muscles. For thinner clients, please use the air cushion head on Level 1 vibration speed to avoid pain or discomfort. For overweight or clients with well-developed muscles, using the round head attachment on Level 2-3 speed improves the treatment effect.

EN

DE


FR

ES

IT

Full Body Massage


ARM & PLAM

The bullet head attachment is recommended for treatments of the palm area, especially on small muscle tissues and joint parts. Levels 1 - 2 is recommended.

Full Body Massage

LEG & FEET

The ball and air cushion attachments are recommended for hip muscles on Levels 2-3 speed to improve the treatment effect. Please avoid using on the spine or tailbone.

EN

DE


FR

ES

IT

Full Body Massage


LEG & FEET


The ball and air cushion head attachments are recommended for thigh muscles on vibration speed Levels 2-3.

Full Body Massage

LEG & FEET

Above mentioned head attachments are also applied for a deeper massage on the back thigh and calf muscles.

The bullet head attachment is recommended for small muscle tissues on the arch and joint parts. Vibration speed level 1 - 2 is therefore recommended.

EN

DE

FR


ES

IT

Full Body Massage

ABDOMEN

During massaging the abdomen, please stay in a standing posture. Level 3 vibration speed is recommended for the abdomen area. The purpose is to use fast shockwaves to vibrate adipose tissue under the skin. Do not use the vibration device on the abdomen continuously for more than 5 minutes to prevent an adverse reaction of the digestive system. Usage on the abdomen is to achieve passive exercise, weight loss, and body shaping treatment. A vibration massage on the abdomen is recommended after sport.


The air cushion head attachment is recommended on the abdomen area and adipose tissue on stomach.

Please do not massage your abdomen immediately after dinner.

Full Body Massage

MORE INFORMATION ABOUT A VIBRATION MASSAGE

No matter which head attachments used, shockwaves of this device are able to vibrate through skin and surface tissues. During the treatment, adding extra pressure on the massage device will increase the intensity and depth of shockwaves therefore improving the treatment without influencing the vibration frequency.

EN

DE


FR

ES

IT

Full Body Massage


MORE INFORMATION ABOUT A VIBRATION MASSAGE

To increase the intensity and depth of shockwaves, it is recommended to hold the muscle parts around the treated area steadily with hands. By preventing the diffusion of shockwaves through skin and surface tissues, the massage gun can achieve a deeper effect on muscle parts.

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